

Customer	Contact	Date	16-Aug-10
----------	---------	------	-----------


Date		Test 2									
Material	pNPP, disodium, hexhydrate	Material condition	.5 mm	Bulk Density (gm/cc)		.1 mm	Bulk Density (gm/cc)				
Desired Sample size	55 mg tgt. wt.	Patriculate size	Crystals >50 micron	Fines	>1%	From Table	Measured	Tap Density	From Table	Measured	Tap Density
Desired accuracy	Range= 50 - +/-10% 60 mg (+/- 5 mg)		rh		53%		0.8	0.85	NA		
Pipette Size	0.125 dia.	Filter Cup:	Filter cup mesh		10 micron						

Micrometer setting	0.40	Sample weight	Control unit settings
	Set	MG	Vacuum (in hg) 12
		54.80	Air (psi) 10
		55.10	
		54.40	
		53.70	
		55.60	
		52.60	
		54.80	
		53.80	
		54.00	
		54.40	
		53.60	
		52.00	
		51.50	
		52.90	
		53.30	
		55.20	
		54.00	
		52.90	
		53.00	
		54.90	
		52.00	
		52.70	

Cycle time (sec)	Aspirate	Level	Dispense	total
	2	0	1	3

53.6381	Av
53.6264	Mean
1.1474	Std Dev
55.6	High
52.0	Low
6.47	%

DISPENSE WEIGHT



Notes/observations

pNPP, disodium, hexhydrate

- 1a. Free flowing white crystals, easily aspirated and dispensed. No issues dispensing into the customers vials.
- 1b. Set up and calibration are straight forward.
- 1c. Requires relatively high air pressure to eject sample from the pipette due to the crystalline shape compacting.
- 1d. Use of the leveling plate is not recommended. Tip pressure comacts crystals impeding the tip discharge. Instead, the tip was tapped to remove excess powder proud of the tip orifice.
- 2.0 Standard 0.125 dia tip used for testing due to small bottle neck dia.
- 3.0 **Aerosol:** No aerosol or loose crystals were observed during or after dispense.
- 4.0 Under normal laboratory conditions, periodic cleaning, based on humidity. Tip/filter cleaning and changeover is <5 minutes.

Summary: Easy material to pipette, even at the small dispense weight. Dispense weight can be normalized by brushing the excess powder proud of the tip orifice