


Micrometer setting
0.10
set

7.9909 Av
7.9908 Mean
0.0406 Std Dev

8.0500 High Range
7.9500 Low
1.2422 %

7.90
8.00
7.97
7.95
7.95
8.00
7.95
8.00
7.95
8.00
7.95
8.00
8.00
8.02
8.00
8.00
8.00

Vacuum (1in hg) 12 Air (psi) 14



Cycle time (sec) Aspirate 1 Level Skip Dispense 1 total 3

7.9824 Av
7.9786 Mean
0.0309 Std Dev

8.02 High Range
7.95 Low
0.87 % See note 3

Notes/observations

1. **Type S:** White, granular, fairly free flowing powder with small percentage of fines <5 micron. High angle of repose.
 - 1a. 8 mg Sample will eject from pipette into syringe at 4 - 6 psi without blowback, aerosol, or flying dust particles.
 - 1a. Dispenses readily into 1.5 ml vial with no blowback aerosol.

2. **PAC:** White, granular powder, uniform particulate size. Low in fines, <50 micron Particulate size. High angle of repose. Tends to agglomerate into a clump in the pipette when aspirated.
 - Had to continually sift the powder to pipette uniform samples.
 - Hydroscopic, Powder gains weight due to moisture if left open .
 - 2.a. Requires fairly high air pressure to fully eject sample from the tip due to accretion in the tip.
 - 2.b. 8.0 gm dispense weight requires 8-9 PSI to fully eject sample into the tube. No aerosol.
 - 2.c. Dispenses readily into 1.5 ml vial with no blowback aerosol.
 - 2.d. Beads carry a large electrostatic charge causing the beads to agglomerate in the supply vessel and pipette tip when dry.
 - 2.e. Recommend static control when dispensing.

3. **PAC test # 2:** Powder was vacuum dried. A 0.062 dia proto tip was fabricated to improve bead ejection during dispense cycle.
 - Hydroscopic, Powder gained weight due to moisture during test. Control sample, 160 mg gained 10 mg (6.3%) @ 30% rh in 30 minutes.
 - Levelling step was eliminated to reduce tip clogging.
 - 3.a. Requires fairly high air pressure to fully eject sample from the tip due to accretion and softening of the beads in the tip.
 - 3.b. 8.0 gm dispense weight requires 12 PSI to fully eject sample into the tube. No aerosol.
 - 3.c. Dispenses readily into 1.5 ml vial with no blowback aerosol.
 - 3.d. Repeatability improves with controlling moisture and use of smaller ID tip.

