


Customer						Date	4/23/02
Material	Canola seed			Est. Bulk Density	0.6		
Desired sample size	100 mg						
Desired accuracy	Plus/Minus 5 mg			Pipette Size	.250 dia	rh	50%
		Sample weight		Control unit settings			
Micrometer setting		gm		Vacuum (in hg)	Air (psi)		
2.50		104.3		5	4		
↓		94.2					
		93.1				wet	
		103.0					
		96.7				wet	
		97.0					
		107.4				wet	
		108.0					
		93.5				wet	
		106.0					
		98.6					
		96.5					
		94.3					
		104.5					
		101.2				refill tray	
		101.1				wet	
		105.3				↓	
		100.3					
		100.3	Av				
		100.2	Mean				
		4.983	Std Dev				



Sample material contained much moisture  
 Small sample size, dried unevenly.

### Sample Size

